

AD-A192 710

SCANNING TUNNELING MICROSCOPY AS A SURFACE CHEMICAL
PROBE(U) MARYLAND UNIV COLLEGE PARK DEPT OF PHYSICS AND
ASTRONOMY E D WILLIAMS 31 MAR 88 AFOSR-TR-88-0427

1/1


UNCLASSIFIED

AFOSR-86-0235

F/G 7/4

NL





②

AD-A192 710

DOCUMENTATION PAGE

Form Approved
OMB No. 0704-0188

1b. RESTRICTIVE MARKINGS DTIC FILE COPY	
2a. SECURITY CLASSIFICATION AUTHORITY	
2b. DECLASSIFICATION/DOWNGRADING SCHEDULE	
4. PERFORMING ORGANIZATION REPORT NUMBER(S)	
5. MONITORING ORGANIZATION REPORT NUMBER(S) AFOSR-TR- 88-0427	
6a. NAME OF PERFORMING ORGANIZATION University of Maryland	6b. OFFICE SYMBOL (if applicable) NC
7a. NAME OF MONITORING ORGANIZATION AFOSR/NC	
7b. ADDRESS (City, State, and ZIP Code) Building 410 Bolling AFB, DC 20332-6448	
8a. NAME OF FUNDING/SPONSORING ORGANIZATION AFOSR	8b. OFFICE SYMBOL (if applicable) NC
9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER AFOSR-86-0235	
10. SOURCE OF FUNDING NUMBERS	
PROGRAM ELEMENT NO. 61102F	PROJECT NO. 2917
TASK NO. A2	WORK UNIT ACCESSION NO.
11. TITLE (Include Security Classification) Scanning Tunneling Microscopy as a Surface Chemical Probe ,	
12. PERSONAL AUTHOR(S) Ellen D. Williams	
13a. TYPE OF REPORT Final	13b. TIME COVERED FROM 86/08/01 TO 88/01/31
14. DATE OF REPORT (Year, Month, Day) 88/03/31	
15. PAGE COUNT 2	
16. SUPPLEMENTARY NOTATION	
17. COSATI CODES	
FIELD	GROUP
18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)	
19. ABSTRACT (Continue on reverse if necessary and identify by block number) <p>This grant sponsored purchase of surface characterization spectrometers, data acquisition hardware, and sample manipulation for an ultrahigh vacuum chamber equipped with a scanning tunneling microscope. This apparatus is designed to map atomic-scale surface morphology of well-controlled samples, particularly for adsorbate-covered metal surfaces. The equipment assembly has been completed and experiments are underway.</p>	
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT <input checked="" type="checkbox"/> UNCLASSIFIED/UNLIMITED <input type="checkbox"/> SAME AS RPT. <input checked="" type="checkbox"/> DTIC USERS	
21. ABSTRACT SECURITY CLASSIFICATION Unclassified	
22a. NAME OF RESPONSIBLE INDIVIDUAL Lt Col Larry W. Buregraf	22b. TELEPHONE (Include Area Code) (202) 767-4960
22c. OFFICE SYMBOL NC	

DTIC ELECTED
MAY 02 1988
Cu E

AFOSR-TR. 88-0427

Final Report

Grant # AFOSR-86-0235

Principal Investigator: Ellen D. Williams

Department of Physics and Astronomy

University of Maryland

College Park, Maryland 20742

March 31, 1988

Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	



88 5_02 134

The instrumentation provided under this grant is in use in our on-going AFOSR research project; "Scanning Tunneling Microscopy as a Surface Chemical Probe". The application of each piece of equipment is listed below.

1. Micro-computer

The micro-computer is dedicated to real-time control and data-acquisition of the STM. Software has been developed to drive the STM scans under computer control while simultaneously reading the output of the feedback circuit. The data are then displayed on the monitor (in low resolution) and stored in an array in memory.

2. Graphics Software

A sophisticated commercial graphics software allows data manipulation and three-dimensional displays after data acquisition is complete.

3. Quadrupole Mass Spectrometer

The mass spectrometer is mounted on the UHV system housing the STM. It is used for partial pressure analysis of the background. In the future it will be used for measurement of thermal desorption spectra following STM scans of adsorbed overlayers.

4. X-Y recorder:

The recorder is used for routine acquisition of Auger and Appearance Potential spectra of samples.

5. LEED system:

The rear-view LEED optics are mounted on the vacuum system housing the STM. They are used for qualitative display of the diffraction structure of the samples. In the future, an existing video detector will be used to acquire quantitative beam profile measurements for comparison with the STM images.

6. Manipulator

The liquid-nitrogen cooled manipulator has been ordered and delivery is scheduled within a month. The cooling feature will be used to monitor the kinetic processes of surface faceting over temperature ranges down to 100 K.

END

DATE

FILMED

6-1988

DTIC